| CANDIDATE | UNIVERSITY OF CAMBRIDGE INTERNATIONAL EXAMINATIONS General Certificate of Education Advanced Subsidiary Level and Advanced Level | ANN. PabaCambride |
|--------------------------|--|--------------------------|
| NAME CENTRE NUMBER | CANDIDATE NUMBER | |
| BIOLOGY | | 9700/34 |
| Advanced Prac | tical Skills 2 | May/June 2012 2 hours |
| Candidates and | wer on the Question Paper. | |
| Additional Mate | rials: As listed in the Confidential Instructions. | |

READ THESE INSTRUCTIONS FIRST

Write your Centre number, candidate number and name on all the work you hand in. Write in dark blue or black ink. You may use a pencil for any diagrams, graphs or rough working. Do **not** use red ink, staples, paper clips, highlighters, glue or correction fluid. DO **NOT** WRITE IN ANY BARCODES.

Answer **all** questions.

You may lose marks if you do not show your working or if you do not use appropriate units.

At the end of the examination, fasten all your work securely together. The number of marks is given in brackets [] at the end of each question or part question.

| For Examiner's Use | | |
|--------------------|--|--|
| 1 | | |
| 2 | | |
| Total | | |

This document consists of **9** printed pages and **3** blank pages.



You are reminded that you have only one hour for each question in the practical examination

You should:

- Read carefully through the whole of Question 1 and Question 2.
- www.papaCambridge.com Plan your use of the time to make sure that you finish all the work that you would like to do.

You will gain marks for recording your results according to the instructions.

1 Plant cells contain an enzyme, catalase, which catalyses the hydrolysis (breakdown) of hydrogen peroxide into oxygen and water. An extract of plant tissue contains catalase.

You are required to:

- make different concentrations of plant extract containing catalase •
- investigate the effect of different concentrations of plant extract by finding the time taken for a square of filter paper, soaked in the plant extract, to rise in hydrogen peroxide solution.

You are provided with:

| labelled | contents | hazard | volume /cm ³ |
|----------|-----------------------------|---------------------|----------------------------|
| Η | hydrogen peroxide solution | harmful irritant | 100 |
| W | distilled water | none | 80 |
| Ρ | 100% plant extract solution | none | 50 |

(a) (i) Decide on the concentrations of plant extract solution you will use in your investigation.

You will need to make up 10 cm³ of each concentration of plant extract solution.

Prepare the space on page 3 to show:

- the concentration of P
- the volumes of **P**
- the volumes of **W**.

Use this space:

You are advised to read steps 1 to 10 before proceeding.

Proceed as follows:

- 1. Prepare the concentrations of plant extract solution as stated in (a)(i).
- 2. Put **H** into the test-tube, filling to within 2 cm from the top.
- 3. Cut squares of filter paper, 1 cm × 1 cm.
- 4. Use forceps to pick up one square of filter paper and dip the whole square into one of the concentrations of plant extract solution in its container.
- 5. Wipe the square against the inside of the container to remove excess plant extract solution from both sides of the square.
- 6. Hold the square so that the top of the square is level with the surface of **H** as shown in Fig. 1.1.



Fig. 1.1

- 7. Release the square (you may need to shake the forceps) and start timing.
- Record the time taken for the square to return to the surface.
 If the time is more than three minutes, stop timing and record 'more than 180'.
- 9. Remove the square from the test-tube.

3

www.papacambridge.com

(ii) Prepare the space below and record your results.

[6]

www.papacambridge.com

(iii) Identify three significant sources of error in your investigation.

| | [3] |
|------|---|
| | |
| (iv) | Suggest how you would make three improvements to this investigation. |
| (iv) | Suggest how you would make three improvements to this investigation. |
| (iv) | Suggest how you would make three improvements to this investigation. |
| (iv) | Suggest how you would make three improvements to this investigation. |
| (iv) | Suggest how you would make three improvements to this investigation. |
| (iv) | Suggest how you would make three improvements to this investigation. |

Catalase catalyses the hydrolysis of hydrogen peroxide into oxygen and water.

A student investigated the effect of changing the concentration of hydrogen peroxide solution on this hydrolysis.

www.papacambridge.com For each concentration of hydrogen peroxide solution the time to collect 25 cm³ of oxygen was recorded.

Table 1.1 shows the results of the student's investigation.

| Та | b | e | 1 | .' | 1 |
|----|---|---|---|----|---|
| | | | | | |

| percentage concentration of hydrogen peroxide | time to collect 25 cm ³ of oxygen /s |
|--|--|
| 4 | 47 |
| 8 | 18 |
| 12 | 15 |
| 16 | 13 |
| 25 | 11 |

(b) (i) Plot a graph of the data shown in Table 1.1.



[4]

| 6 Mana Baba | |
|--|---|
| Explain the effect of changing the concentration of hydrogen peroxide. | |
| | oria |
| | .9e.C. |
| | OT |
| | |
| | |
| | |
| | |
| [3] | |
| [Total: 22] | |
| | 6 Explain the effect of changing the concentration of hydrogen peroxide. |

2 K1 is a slide showing transverse sections through blood vessels.

Select one large artery and one large vein.

(a) (i) Draw a large plan diagram of two different types of blood vessel on slide K1.Label the muscle layer in the artery.

[5]

7

Fig. 2.1 is a photomicrograph of a transverse section through lung tissue.





(ii) Make a large drawing to show the walls of **two** whole alveoli that are touching each other.

Label the gas exchange membrane.

www.papacambridge.com

Fig. 2.2 is a photomicrograph of a transverse section through healthy lung tissue.

Fig. 2.3 is a photomicrograph of a transverse section through unhealthy lung tissue.





magnification ×95

Fig. 2.2

magnification ×95

Fig. 2.3

To find the mean diameter of an alveolus in the healthy lung, a student measured five alveoli, as shown by the lines C, D, E, F and G in Fig. 2.2.

Use the magnification to calculate the mean actual diameter, in μ m, of an alveolus (iii) using the lines C, D, E, F and G.

You may lose marks if you do not show your working or if you do not use appropriate units.

.....μm [3]

8

- 9
- www.papacambridge.com (iv) Prepare the space below so that it is suitable for you to show the obs differences between the specimens in Fig. 2.2 and Fig. 2.3.

Record your observations in the space you have prepared.

[5]

[Total: 18]



BLANK PAGE

10



BLANK PAGE

11



BLANK PAGE

Copyright Acknowledgements:

 Question 2, Fig. 2.1
 ED RESCHKE/PETER ARNOLD INC./SCIENCE PHOTO LIBRARY

 Question 2, Fig. 2.2
 BIOPHOTO ASSOCIATES/SCIENCE PHOTO LIBRARY

 Question 2, Fig. 2.3
 BIOPHOTO ASSOCIATES/SCIENCE PHOTO LIBRARY

Permission to reproduce items where third-party owned material protected by copyright is included has been sought and cleared where possible. Every reasonable effort has been made by the publisher (UCLES) to trace copyright holders, but if any items requiring clearance have unwittingly been included, the publisher will be pleased to make amends at the earliest possible opportunity.

University of Cambridge International Examinations is part of the Cambridge Assessment Group. Cambridge Assessment is the brand name of University of