

**CANDIDATE** 

### UNIVERSITY OF CAMBRIDGE INTERNATIONAL EXAMINATIONS General Certificate of Education Advanced Subsidiary Level and Advanced Level

NAME				
CENTRE NUMBER		CANDIDATE NUMBER		
CHEMISTRY Advanced Pract	ical Skills 2	Octo	9701/36 ber/November 2012	
Candidates ans	wer on the Question Paper.		2 hours	

# **READ THESE INSTRUCTIONS FIRST**

Write your Centre number, candidate number and name on all the work you hand in. Give details of the practical session and laboratory where appropriate, in the boxes provided. Write in dark blue or black pen.

As listed in the Confidential Instructions

You may use a soft pencil for any diagrams, graphs or rough working. Do not use staples, paper clips, highlighters, glue or correction fluid. DO NOT WRITE IN ANY BARCODES.

Answer all questions.

Additional Materials:

You may lose marks if you do not show your working or if you do not use appropriate units. Use of a Data Booklet is unnecessary.

Qualitative Analysis Notes are printed on pages 15 and 16.

At the end of the examination, fasten all your work securely together. The number of marks is given in brackets [ ] at the end of each question or part question.

Session				
Laboratory				

For Examiner's Use			
1			
2			
3			
Total			

This document consists of 13 printed pages and 3 blank pages.



**BLANK PAGE** 

www.papaCambridge.com You are to determine the concentration of aqueous copper(II) sulfate by titration concentration of Cu2+ ions in a solution can be found by reaction with an excess of aquit iodide ions to produce iodine. The amount of iodine formed can be found by titration w thiosulfate ions,  $S_2O_3^{2-}$ .

**FB 1** is aqueous copper(II) sulfate, CuSO<sub>4</sub>.

FB 2 is 0.100 mol dm<sup>-3</sup> sodium thiosulfate, Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>.

FB 3 is aqueous potassium iodide, KI. starch indicator

Read through the instructions carefully before starting any practical work.

#### (a) Method

- Fill the burette with **FB 2**.
- Pipette 25.0 cm<sup>3</sup> of **FB 1** into a conical flask.
- Use a measuring cylinder to add 10 cm<sup>3</sup> of **FB 3** into the conical flask.
- Titrate this mixture with FB 2 until the colour of the mixture changes from brown to yellow-brown. An off-white precipitate will also be present in the flask throughout the titration.
- Add approximately 1 cm<sup>3</sup> of starch indicator.
- Continue the titration until the blue-black colour of the starch-iodine complex just disappears leaving the off-white precipitate.
- Perform a rough titration and record your burette readings in the space below.

The rough titre is	3	cm <sup>3</sup>
--------------------	---	-----------------

- Carry out as many accurate titrations as you think necessary to obtain consistent results.
- Make certain any recorded results show the precision of your practical work.
- Record in a suitable form below all of your burette readings and the volume of FB 2 added in each accurate titration.

I	
II	
III	
IV	
V	
VI	

[6]

(b) From your accurate titration results, obtain a suitable value to be used in your calculate. Show clearly how you obtained this value.

www.PapaCambridge.com 25.0 cm<sup>3</sup> of **FB 1** required ...... cm<sup>3</sup> of **FB 2** [1]

## (c) Calculations

The equations for the formation of iodine and its reaction with thiosulfate ions are given below.

$$2Cu^{2+} + 4I^{-} \rightarrow 2CuI + I_{2}$$

$$I_2 + 2S_2O_3^{2-} \rightarrow S_4O_6^{2-} + 2I^{-}$$

Show your working and appropriate significant figures in the final answer to each step of your calculations.

(i) Calculate the number of moles of thiosulfate ions,  $S_2O_3^{2-}$ , present in the volume of FB 2 in (b).

moles of  $S_2O_3^{2-} = \dots$  mol

(ii) Using the equations above, deduce the number of moles of Cu2+ ions present in each 25.0 cm<sup>3</sup>.

moles of  $Cu^{2+} = \dots mol$ 

(iii) Calculate the concentration, in mol dm<sup>-3</sup>, of copper(II) sulfate in **FB 1**.

concentration of  $CuSO_4 = \dots mol dm^{-3}$ 

[3]

(d) Three students repeated the experiment but each obtained different values concentration of CuSO<sub>4</sub>.

The students each suggested possible improvements.

www.PapaCambridge.com Student 1 suggested that a larger quantity of starch should be added. Student 2 suggested that a larger volume of potassium iodide, **FB 3**, should be added. Student 3 suggested that the contents of the conical flask should be filtered before

titration.

Comment on the effectiveness of each of these possible improvements. Explain your answers.

Student 1			
Student 2			
Student 3			
	 	 	[2]

[Total: 12]

www.PapaCambridge.com You are to determine the enthalpy change for the reaction between aqueous con 2 sulfate and zinc. The enthalpy change of reaction can be found by measuring the tempera change when powdered zinc is added to aqueous copper(II) sulfate.

FB 4 is 1.10 mol dm<sup>-3</sup> aqueous copper(II) sulfate, CuSO<sub>4</sub>. powdered zinc

#### (a) Method

- Weigh a 100 cm<sup>3</sup> beaker.
- In the beaker weigh out between 2.1 g and 2.3 g of powdered zinc.
- Record the weighings and the mass of zinc in the space below.

mass of zinc used = ..... g

- Support the plastic cup in a 250 cm<sup>3</sup> beaker.
- Use a measuring cylinder to transfer 50 cm<sup>3</sup> of **FB 4** into the plastic cup.
- Measure and record in the table below, the initial temperature of **FB 4** in the cup.
- Start the stop watch. Measure and record the temperature of FB 4 in the cup after 1 minute, 2 minutes and 3 minutes.
- At time 3½ minutes, add the weighed zinc to **FB 4** in the cup and stir the mixture.
- From time 4 minutes, continue to stir the mixture and measure the temperature of the contents of the cup to complete the table.

#### Results

time/min	0	1	2	3	4	5	6	7	8	9	10	11	12
temperature/°C													

[2]

(b) (i) On the axes opposite, plot the temperature (y-axis) against time (x-axis). The temperature axis should allow you to include a point at least 5 °C greater than the maximum temperature recorded.

www.PapaCambridge.com

(ii) Complete the graph to show how the temperature of the contents of the cup varies with time.

- Draw one straight line through the points between time 0 minutes and 3 minutes.
- Draw one straight line through the points after the maximum was reached.
- Extrapolate these two lines and draw a vertical line at time 3½ minutes.

#### (c) Calculation

(i) Use your graph to determine the change in temperature at 3½ minutes.

www.PapaCambridge.com change in temperature = .....°C

(ii) Calculate the heat energy produced in the reaction. (You may assume that 4.3 J are required to raise the temperature of 1.0 cm<sup>3</sup> of any solution by 1.0 °C.)

heat energy produced = ...... J

(iii) The reaction between zinc and aqueous copper(II) sulfate is a displacement reaction shown in the equation below.

$$Zn(s) + Cu^{2+}(aq) \rightarrow Cu(s) + Zn^{2+}(aq)$$

From the mass of zinc added and the concentration of **FB 4**, show that the copper(II) sulfate was in excess in your reaction.

[A,: Zn, 65.4]

(iv) Assuming that the copper(II) sulfate was in excess, use your answer to (ii) to calculate the enthalpy change of the reaction between Zn(s) and Cu<sup>2+</sup>(aq).

Give you answer in kJ mol<sup>-1</sup> and include the relevant sign.

enthalpy change of reaction = ...... ..... kJ mol<sup>-1</sup> sign value

[6]

. /	-or	
•	niner's	
	ę	

(d) One source of error in this experiment is due to the accuracy to which the therm can be read.

www.PapaCambridge.com What is the maximum error in a single temperature reading on a thermometer with graduations at 1°C?

maximum error = ..... °C

Calculate the maximum percentage error when measuring a temperature rise of 12.0 °C.

maximum percentage error = ...... % [2]

[Total: 14]

At each stage of any test you are to record details of the following.

- colour changes seen
- the formation of any precipitate
- the solubility of such precipitates in an excess of the reagent added

Where gases are released they should be identified by a test, **described in the appropriate place in your observations**.

You should indicate clearly at what stage in a test a change occurs. Marks are **not** given for chemical equations.

No additional tests for ions present should be attempted.

If any solution is warmed, a boiling tube MUST be used.

Rinse and reuse test-tubes and boiling tubes where possible.

Where reagents are selected for use in a test, the name or correct formula of the element or compound must be given.

You are provided with solutions FB 5, FB 6, FB 7 and FB 8. FB 5 and FB 6 each contain a compound of a transition element.

Half fill a 250 cm³ beaker with water. Heat to approximately 80 °C, then stop heating and switch off the Bunsen burner. You will need this as a hot water bath on **(b)(i)**. Continue work on **(a)** while the water heats.

#### (a) (i) Carry out the following tests on FB 5.

test	observations
To 1 cm depth of <b>FB 5</b> in a test-tube, add aqueous ammonia.	
To 1 cm depth of <b>FB 5</b> in a test-tube, add aqueous sodium hydroxide.	
To 1 cm depth of <b>FB 5</b> in a test-tube, add aqueous barium chloride or aqueous barium nitrate then,	
add an excess of either hydrochloric acid or nitric acid.	

For miner's e

(1)	Carry out the following tests on <b>FE</b>	s, if any, can you reach about the identity of F
	test	observations
boil	1 cm depth of <b>FB 6</b> in a ing tube, add 1 cm depth <b>B 7</b> then,	
Plac war	1 cm depth of ethanol. ce the boiling tube in the m water bath and leave for w minutes.	
	1 cm depth of <b>FB 6</b> in a -tube, add 1 cm depth of <b>8</b> .	
(ii) [	From these tests suggest identities	s for the following
	The anion in <b>FB 6</b> is	
	The cation in <b>FB 7</b> is	
٦	The cation in <b>FB 8</b> could be	or
	Suggest a test to determine which Do not carry out this test.	of the two possible cations is present in <b>FB 8</b> .
		[-7]
•		[7]

	Q.	
(d)	Suggest what happened to the ethanol when it was warmed with the mixture of <b>F FB 7</b> .	For miner's e
(e)	You are to devise and carry out a test to confirm the identity of the cation in <b>FB 7</b> .	Se. COM
	Record the test you use and the results of the test in the space below.	

[1]

[Total: 14]

**BLANK PAGE** 

**BLANK PAGE** 

## **Qualitative Analysis Notes**

Key: [ppt. = precipitate]

## 1 Reactions of aqueous cations

ia	reaction with						
ion	NaOH(aq)	NH <sub>3</sub> (aq)					
aluminium, Al³+(aq)	white ppt. soluble in excess	white ppt. insoluble in excess					
ammonium, NH <sub>4</sub> <sup>+</sup> (aq)	no ppt. ammonia produced on heating	_					
barium, Ba <sup>2+</sup> (aq)	no ppt. (if reagents are pure)	no ppt.					
calcium, Ca <sup>2+</sup> (aq)	white ppt. with high [Ca <sup>2+</sup> (aq)]	no ppt.					
chromium(III), Cr³+(aq)	grey-green ppt. soluble in excess giving dark green solution	grey-green ppt. insoluble in excess					
copper(II), Cu <sup>2+</sup> (aq)	pale blue ppt. insoluble in excess	blue ppt. soluble in excess giving dark blue solution					
iron(II), Fe <sup>2+</sup> (aq)	green ppt. turning brown on contact with air insoluble in excess	green ppt. turning brown on contact with air insoluble in excess					
iron(III), Fe³+(aq)	red-brown ppt. insoluble in excess	red-brown ppt. insoluble in excess					
lead(II), Pb <sup>2+</sup> (aq)	white ppt. soluble in excess	white ppt. insoluble in excess					
magnesium, Mg <sup>2+</sup> (aq)	white ppt. insoluble in excess	white ppt. insoluble in excess					
manganese(II), Mn²+(aq)	off-white ppt. rapidly turning brown on contact with air insoluble in excess	off-white ppt. rapidly turning brown on contact with air insoluble in excess					
zinc, Zn <sup>2+</sup> (aq)	white ppt. soluble in excess	white ppt. soluble in excess					

[Lead(II) ions can be distinguished from aluminium ions by the insolubility of lead(II) chloride.]

#### 2 Reactions of anions

ion	reaction
carbonate, CO <sub>3</sub> <sup>2-</sup>	CO <sub>2</sub> liberated by dilute acids
chromate(VI), CrO <sub>4</sub> <sup>2-</sup> (aq)	yellow solution turns orange with H <sup>+</sup> (aq); gives yellow ppt. with Ba <sup>2+</sup> (aq); gives bright yellow ppt. with Pb <sup>2+</sup> (aq)
chloride, C <i>l</i> <sup>-</sup> (aq)	gives white ppt. with Ag <sup>+</sup> (aq) (soluble in NH <sub>3</sub> (aq)); gives white ppt. with Pb <sup>2+</sup> (aq)
bromide, Br <sup>-</sup> (aq)	gives cream ppt. with Ag <sup>+</sup> (aq) (partially soluble in NH <sub>3</sub> (aq)); gives white ppt. with Pb <sup>2+</sup> (aq)
iodide, I <sup>-</sup> (aq)	gives yellow ppt. with Ag <sup>+</sup> (aq) (insoluble in NH <sub>3</sub> (aq)); gives yellow ppt. with Pb <sup>2+</sup> (aq)
nitrate, NO <sub>3</sub> -(aq)	NH <sub>3</sub> liberated on heating with OH <sup>-</sup> (aq) and A <i>l</i> foil
nitrite, NO <sub>2</sub> -(aq)	$NH_3$ liberated on heating with $OH^-(aq)$ and $Al$ foil; NO liberated by dilute acids (colourless $NO \rightarrow$ (pale) brown $NO_2$ in air)
sulfate, SO <sub>4</sub> <sup>2-</sup> (aq)	gives white ppt. with Ba <sup>2+</sup> (aq) or with Pb <sup>2+</sup> (aq) (insoluble in excess dilute strong acids)
sulfite, SO <sub>3</sub> <sup>2-</sup> (aq)	SO <sub>2</sub> liberated with dilute acids; gives white ppt. with Ba <sup>2+</sup> (aq) (soluble in excess dilute strong acids)

#### 3 Tests for gases

gas	test and test result
ammonia, NH <sub>3</sub>	turns damp red litmus paper blue
carbon dioxide, CO <sub>2</sub>	gives a white ppt. with limewater (ppt. dissolves with excess CO <sub>2</sub> )
chlorine, Cl <sub>2</sub>	bleaches damp litmus paper
hydrogen, H <sub>2</sub>	"pops" with a lighted splint
oxygen, O <sub>2</sub>	relights a glowing splint
sulfur dioxide, SO <sub>2</sub>	turns acidified aqueous potassium dichromate(VI) from orange to green

Permission to reproduce items where third-party owned material protected by copyright is included has been sought and cleared where possible. Every reasonable effort has been made by the publisher (UCLES) to trace copyright holders, but if any items requiring clearance have unwittingly been included the publisher will be pleased to make amends at the earliest possible opportunity.

University of Cambridge International Examinations is part of the Cambridge Assessment Group. Cambridge Assessment is the brand name of University of